## English

### ST-3 Sensor TESTER, Alice

**Intended Use**

The ST-3 Alice is intended to test sensors, leads, and microphones used in sleep studies and related areas.

**Package Contents**

- ST-3 Alice Sensor Tester Module
- Continuity Test Cable
- 9 Volt Alkaline Battery
- User's Guide

Square Wave output cables must be ordered separately.

**Caution**

This product is for equipment testing purposes only and is not to be used as a diagnostic patient device, apnea monitor, or in a life supporting, or life sustaining situation.

**Warning**

Sensors and equipment must be disconnected from patient before testing. This device cannot yield impedance testing of sensors that have been applied to patients. At no time should this device be connected to any line power (115-120 Volts).

**WEEE-RoHS/Disposal**

⚠ If applicable, separate collection for electrical and electronic equipment per EC Directive 2002/96/EC. Otherwise, dispose of this device in accordance with local regulations.

| Specifications   |  |
|--|--|
| Size   | 70mm(W) x 126mm(L) x 24mm(H)   |
| Weight (with battery), ST-3 Alice  | 130 grams  |
| Power  | Battery, 9 volt alkaline   |
| Square Wave Output Amplitude & Frequency   | 100 microvolts peak-to-peak, polarity changes every 2 seconds (~0.25Hz) simulating EEG machine calibration outputs. When using high level DC output cables this will be 0 – 1 volt output level at the same frequency. |
| Storage Temperature  | -20°C (-4°F) to 60°C (140°F)   |
| Operating Temperature  | 5°C (41°F) to 40°C (104°F)   |
| Operating and Storage Humidity   | 15 - 95% Non-condensing  |
| <b>Output Cable Options</b>  |  |
| The square wave output generator may be adapted to different inputs for the various devices in a sleep testing facility. To do this there are cables with different plugs available for these connections: |  |
| Part Number  | Connection Type  |
| 1075   | 1.5mm safety connectors  |
| 1076   | 1/8" male phone plug   |
| 1077   | 3/32" male phone plug  |
| 1079   | 1/8" female phone receptacle   |
| 1081   | 1/4" male phone plug   |

**Calibration and Maintenance**

**Calibration**

The ST-3 Alice requires no calibration. It is designed to be a qualitative device.

**Maintenance**

There are no operator-serviceable parts in the ST-3 Alice except the 9 volt battery. Inspection of battery and cable connectors is useful in maintaining proper operation of the ST-3 Alice.


Illustration 1

**Inputs and Display**

**Selector Switch (1)**

This is a 4 position switch:

Other is for testing all standard sensors with 1.5mm safety connectors Mic/Therm is for testing microphones and thermistors Effort is for testing all piezo effort sensors Off turns power to the unit off and preserves battery life

**Alice 3 Effort (2)**

This plastic lemo connector is for testing Alice 3 effort sensors

**Alice 3 Mic/Therm (3)**

This 3/32" (2.5mm) phono input is for testing thermistors and microphones for the Alice 3 system

**Alice 4 (4)**

This triaxial lemo connection is for testing effort sensors, microphones, and thermistors made for the Alice 4 system

**Alice 5 (5)**

To test Alice 5 effort sensors, the switch must be in effort position with sensor in Alice 5 inputs. To test Alice 5 thermistors, the switch must be in Mic/Therm position

**Sensor Test (5)**

The two 1.5 mm safety receptacle inputs on the face of the tester are for testing all standard sensors (with selector switch in Other position)

**Continuity Test (6)**

Plug in the cable which terminates with a 1.5mm safety receptacle and alligator clip. This cable is used to test single wire electrodes to make certain there are no breaks in the wire.

**Square Wave (7)**

Output cables (not included) are plugged into this port to inject the square wave output signal into PSG data acquisition units to test cable connections and jackboxes. The type of output cable is determined by the inputs necessary for interfacing with the specific data acquisition system in use. \*Set the low

frequency filter of the system to the lowest position (time constant to its highest position) to observe the square wave at its best resolution. Amplitude will be +/- 50 microvolts; with high level DC output cables this will be 0 – 1 volt output.

**Display (8)**

The display is a LED bar graph. Sensor excitation determines the movement of the bar. As an effort sensor is stretched the light will travel back and forth, as a snore sensor is tapped it will bounce quickly, as a thermocouple is heated it will move in one direction and as it cools it will move in the other.

**Removing The Battery**


Gently slide the battery cover down and away from the tester as indicated on the cover. Hold the tester in one hand. Using the fingers of the other hand, gently push the bottom of the battery toward the right side of the tester and lift it out.

**Inserting A New Battery**

Gently slide the battery cover down and away from the tester as indicated on the cover. Check the positive (+) and negative (-) polarity indicated on the bottom of the battery compartment before inserting the battery. Hold the tester in one hand. With the fingers of the other hand, snap the battery into place against the power contacts. Gently slide the battery cover up.

**Warranty**

Pro-Tech Services, Inc. warrants the ST-3 Sensor Tester module to be free of defects in materials and workmanship for a period of one year from the date of sale by Pro-Tech Services, Inc. If any product fails to perform in accordance with this warranty, the sole liability of Pro-Tech Services, Inc. is limited to repair or, at its option, replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Services, Inc. be liable for any indirect, special or consequential loss or damage, including economic loss, loss of business or overhead, or any punitive damages, arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Services, Inc. cannot be held liable for conditions resulting from damage, marginal performance or malfunctions caused by misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories, or any other defect not related to materials or workmanship. This warranty is in lieu of all other warranties, expressed or implied, and is extended only to the original purchaser. Features and specifications are subject to change without notice.

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